





iHILIC[®]

Advancing HILIC Separations in UHPLC and HPLC



- Zwitterionic charge modulated amide and diol HILIC columns
- Complementary selectivities for separation of polar compounds
- Excellent durability and ultra-low bleeding
- Universal columns for LC-MS based "Omics" studies and other applications
- iHILIC[®]-Fusion and iHILIC[®]-Fusion(+):
1.8, 3.5, and 5 μm; pH 2-8
- iHILIC[®]-Fusion(P) and iHILIC[®]-(P) Classic:
5 μm; pH 1-10

